



Application of Laser Diffraction with Dry Dispersion for Automatic Particle Size Analysis

U. Kesten *

presented during

**3 ° Congresso Brasileiro de Cimento
Sao Paulo, 8. - 12. novembro 1993**

* Dr.-Ing. U. Kesten
Sympatec GmbH
Burgstätter Str. 6
D - 38678 Clausthal

Application of Laser Diffraction with Dry Dispersion for Automatic Particle Size Analysis

U. Kesten

Introduction

The characterisation of the particle size distribution of cement has become one of the important parameters for the determination and the control of quality. With the introduction of the laser diffraction principle into industrial quality control some 15 years ago the basis for further development was laid.

Sympatec in 1984 introduced a new generation of laser diffraction analysers, based on the simple idea that the instruments used for the characterisation of quality parameters have to adapt to the products and processes and not vice versa, i.e. that the products have to adapt to the instruments, which was the common way of how analytical equipment was designed.

This basic rule consequently led to the development of a range of laser diffraction analysers and dispersing systems that are of modular design thus being easily adaptable to a wide variety of different products and processes.

Dry powders have to be dispersed dry, wet products and suspensions have to be dispersed in the wet state. This is the practical and consequent conclusion.

Sympatec GmbH being founded in 1984 in Germany have developed a to date unique instrument for the dry dispersion of powders in the particle size range from below $0.5 \mu\text{m}$ to more than 3 mm. This patented instrument, the RODOS, allows for a smooth, complete and product adaptable dispersion.

In combination with the laser diffraction instrument, HELOS, a powerful system for industrial particle size analysis of the powders is available.

Company Profile and Technical Highlights

The development of Sympatec GmbH can be taken from the company profile (fig. 1). Sympatec, standing for SYstem PArTicle TECnology, since its foundation in 1984 as a hightech spinout of the University of Technology Clausthal, Germany (famous for its mining tradition), has featured remarkable Highlights in the development of the laser diffraction technology (fig. 2). This certainly underlines Sympatec's claim for technical leadership.

Laser Diffraction Principle

The history of Laser Diffraction as an advanced optical analysis principle for fast determination of particle size distributions in a wide size range, i.e. from 0.1 to $x.000 \mu\text{m}$, begins in the early 19th century with the experimental work of Mr. Fraunhofer in Munich (fig. 3).

Those experiments are the basis of what is known today as Fraunhofer diffraction physics.

From these first experiments it took about 130 years for the technical realisation of the first laser diffraction particle size analyser in 1972 [1].

The optical setup for the generation of diffraction patterns consists of a laser as light source, an optical beam expander that widens the very fine laser beam to more than 10 mm in diameter, and a measuring zone where the sample is fed and taken out.

The measuring zone is the position where the interactions between the monochromatic light and the particles take place and the diffraction patterns are generated (fig. 4).

These phenomena of intensity distribution are focussed on the surface of an optical detector with the help of a Fourier transformation lens of a distinctive focal distance. The detector usually consists of up to 31 semicircular elements. They are transforming the optical light intensity into the corresponding electric current (fig. 5).

A mathematical relationship can be established between the diameter of the particles and the intensity patterns [2]. For a particle size distribution, consisting of fine, intermediate and coarse particles, the diffraction pattern monitored with the detector, comprises a two parameter information; it is the distribution of particle sizes and the quantities (fig. 6).

An extremely important feature that determines the quality of the measurement to a very large extent is the accurate alignment of the optical setup, i.e. the precise focussing of the central laser beam to the centre of the detector. Sympatec have realised an focussing system that automatically and constantly checks this alignment and if necessary corrects it before a measurement can be released (fig. 7).

In a cutaway view of the HELOS/RODOS system the important elements of a technical realisation are shown (fig. 8).

Sympatec Assortment

As a result of the idea of adaptation of the instruments to the products and processes a completely modular conception of the Sympatec product range was realised.

The Sympatec assortment comprises instruments for off-, auto- and on-line particle size analysis with laser diffraction (fig. 9).

The HELOS Laser Diffraction Sensor is manufactured in four different types for a size range between 0.1 μm and 3.500 μm .

Dispersing Systems are provided for dry and wet applications in up to more than 20 variations.

Analytical Software allows for the processing of particle size analysis data. Sample couplers are used for the connection of the analysis system to the production lines if a complete on-line particle size analysis for continuous quality control is required.

Independent of the type of HELOS sensor all instruments show the same technical features. The Central Unit, made of high precision mechanical components, houses also the optical and electrical elements. The Systems Computer, being integral part of the instrument, underlines Sympatec policy of using high quality parts for best system performance (fig. 10).

A typical HELOS/RODOS configuration for off-line application including the Process Control System and a printer, is shown in the photo (fig. 11).

The variety of Dispersing Systems for dry and wet dispersion allows for the best choice with regard to the products to be analysed. The Matrix of dispersers shown can be used for the selection of the appropriate system (fig. 12).

The RODOS, standing for Rotating Dosing and Dispersing System, is a patented disperser for dry powders [3]. Due to its unique two stage design it can be used with powders in a very wide size range from below 0.5 μm to more than 3 mm (fig. 13).

Depending on the dispersion characteristics different adaptors are applied to feed the sample to the gas-solids injector, being the centre piece of the disperser.

For very fine and sticky powder the sample is fed in surplus into a rotary groove with the help of a vibratory feeder. Additional scrapers and rollers are compacting the powder in the groove. Thus a homogenous mass flow of powder is presented to the rear end of the disperser. The gas solids injector is driven by compressed air. In a nozzle it is accelerated to velocity of about 50 m/s.

Due to the acceleration of the compressed air, a vacuum is produced at the rear end of the disperser and the previously homogenized powder is taken into the disperser. In a zone of high velocity gradient between the air and the agglomerated particles, shear forces are applied to destroy these particle aggregates. Additional particle-to-wall and particle-particle collisions and centrifugal forces generated by the rotation of such aggregates completely and smoothly take the agglomerates apart, down to the primary particle sizes.

The powder leaves the disperser as an aerosol in a free jet and crosses the laser beam of the HELOS instrument at a very short distance from the exit of the disperser. After passage of the laser beam the aerosol is completely taken out of the air with the help of an extraction nozzle and a vacuum cleaner.

The advantages of the dry dispersion with RODOS over suspension (wet) dispersion are manifold and include complete and proper product specific dispersion, adaptable dispersing forces and variable sample size, highest measurement frequency available to day, as well as a guaranteed long life time of the disperser with minimum operating costs (fig. 14).

With regard to application in the cement industry RODOS is the natural way for the dispersion of

- cement
- raw meal
- coal dust
- fly ash
-

Automation

In addition to the off-line instruments designed for application in labs of all branches with manually operated dispersing systems, Sympatec have created a new line of instruments, called auto-line, for application in automatic laboratories for continuous quality control. These instruments automatically perform a particle size analysis once a sample has been filled into the disperser. They can be operated as stand alone solutions and, due to their ability for cooperation with robots and central process control systems, they can be integrated into the laboratory automation systems supplied by well known plant manufacturers.

A typical cycle time for a measurement with an integrated Sympatec auto-line system with the automatic dry dispersing unit auto-RODOS module totals to less than two minutes (fig. 15).

Included in the cycle is the information exchange between the central control computer of the automatic lab and the HELOS systems computer, the fill in of the sample, the dispersion and measurement of the size distribution, the evaluation of data and transfer to the host computer [4].

The most significant differences in the differential distribution curves are those between the charge 1 and charge 2: Charge two is showing a remarkably higher amount of coarse particles (fig. 23).

With regard to the long time stability of the dry dispersing system auto-RODOS, Sympatec guarantees a life time of 25.000 measurements of PZ 35 F, i.e. a life time of typically one year in an automatic production control laboratory. The standard deviation over the life time is less than 1 % or expressed in terms of the BLAINE surface, inferior to +/- 40 cm²/g for 4000 cm²/g (fig. 24).

Compared to the operating costs of systems that are using suspension dispersion in alcohol, the dry dispersing system RODOS has an advantage of probably more than 50%.

Conclusion

As a conclusion the following important facts can be stated.

HELOS and auto-RODOS module are the superior components for automatic production control of dry powders they

- * represent the state of the art of laser diffraction technology
- * have been established world wide over the last years
- * are using a unique dry dispersing system
- * allow for product specific dispersion of powders
- * are designed for a high measurement frequency of up to 200 measurements/day
- * have a guaranteed life time of 25000 measurements
- * are thus allowing for lowest operating costs.

References

- [1] *Stephan Röthele, Ulrich Kesten:*
The Application of Laser Diffraction Technology,
2nd World Congress of Particle Technology 1990,
19. - 22. September 1990, Kyoto, Japan
- [2] *Michael Heuer, Kurt Leschonski:*
Results Obtained with a New Instrument for the Measurement
of Particle Size Distributions from Diffraction Patterns,
Part. Charact. 2 (1985) 7 - 13
- [3] *Kurt Leschonski, Stephan Röthele, Ulrich Menzel:*
A Special Feeder for Diffraction Pattern Analysis
of Dry Powders,
Part. Charact. 1 (1984) 161 - 166
- [4] *Stephan Röthele, Heinz Naumann, Udo Brandis:*
on-line Sampling and Sample Splitting -
Principle and Instrumentation
6th Conference Granulometry, 21. - 22. Nov. 1989, Dresden



Company profile Sympatec GmbH SYsteM-PARticle-TEChnology

- ★ Foundation 1984
- ★ Location in Germany Clausthal (Sales, Administration, After-Sales)
Goslar (R & D, Distribution, Manufacturing)
- ★ Venture Capital Swiss Holding Company
- ★ Products & know-how Powder Technology Department
TU Clausthal, University of Karlsruhe (TH)
- ★ Awards
 - 1985: technology transfer award of
IHK Braunschweig
(chamber of commerce);
dry dispersion technology: RODOS
 - 1986: technology transfer award of the
German Federal Minister of
Education and Science;
on-line-particle-size-analysis
 - 1992: Arnold-Euken-Award;
GVC: **OPUS ULTRASONIC
EXTINCTION**
- ★ Marketing strategy
 - laser diffraction for dry
dispersed bulk materials
to below 1 μm
 - development of dispersing systems
for the adaptation of laser diffraction
to products
 - modular conception of instruments
- ★ Manufacturing strategy
 - lean production
 - assembly-plant
 - high-tech-product
- ★ Employees 48 in 1993

Highlights of technical leadership

- 1984 RODOS adaptation of laser diffraction for dry analysis from
less than 1 μm , with patented two-stage dispersion
system
- HELOS highly accurate for auto-aligned 31 classes, combined
with parameter-free Philips-Twomey solution
- 1985 QX programme package for universal processing of particle
size distributions (psd)
- 1986 ROPRON patented sample-coupler with two-stage in-stream sample
splitter (dry processes)
- 1987 SUBMICRON pioneering application of Fraunhofer diffraction
from 0.1 μm , with 11 points measured below 1 μm
- 1988 REMO/PARA Software modules for comparison of different psa-methods,
for example, laser diffraction with sieve analysis
- 1989 auto-RODOS automatic versions of dispersing systems, computer
controlled, roboter-adaptable
- 1990 auto-SUCELL
auto-GRADIS
- TRIMO dynamic trigger-ignition
- SAFIR Sampling Finger Robot (wet processes)
- 1991 QT Quality time-accelerator
- Spray-Sizer universal adaptor for spray-applications
- MAGIC HELOS-Sensor with auto-ranging
- OPUS acoustic on-line-psa using ultra sonic-extinction
- 1992 auto-RODOS-module integrated dry disperser with extremely extended
long-term stability
- TOPMICRON range extension with deflection amplifier
- MIE high speed tool with unlimited refractive index stock


Laser Diffraction

History


Advanced analysis principle for fast determination of particle size distributions in a wide size range i.e. from 0,1 - X.000 μm

Elements	Definition	Year
Theory	Fraunhofer diffraction physics	1840
Method	He-Ne-laser as high energy, coherent light for the generation of particle-light interaction	1960
Sensor	multi-element semiconductor first used as aiming device for military application	1965
Processing	high performance microcomputers	1970

Result	particle size analysis with laser diffraction	since 1972
--------	------------------------------------------------------	------------




optical set up for the generation of diffraction patterns



**Multi-element detector
with 31 classes**

3E08K10U3E



**Intensity distribution of
Fraunhofer diffraction patterns
in the focal plane**



**single particle of
diameter x**

**particle size
distribution $q_0(x)$**

3E08K10U2E



cutaway view of
HELOS/RODOS-system

- 1 He-Ne laser with beam expander
- 2 optical bench system
- 3 dry dispersing system RODOS
- 4 deflection mirrors
- 5 measuring volume
- 6 Fourier lens
- 7 multi-element detector with autofocus-system



SYMPATEC


instruments for off-, auto- and on-line particle size analysis

- ★ **Laser Diffraction : Sensors** : for the analysis of particle size distributions
 - HELOS/BA* : 0.1 - 350 μm
 - KA* : 0.1 - 1750 μm
 - LA* : 0.1 - 3500 μm
 - VARIO* : 0.5 - 3500 μm
 - (* also with MAGIC option)

- ★ **Dispersing Systems** : for the adaptation of the sensors to powders, suspensions, emulsions, sprays
 - dry: **RODOS***, **GRADIS***
 - wet: **SUCELL***, **CUVETTE***
 - (* also in auto-line version)

- ★ **Analytical Software** : for processing of psa data and correlation with other analysis principles
 - QX, REMO & PARA,**
 - TRIMO, QT**

- ★ **Sample Couplers** : for the connection of the analysis systems to the production lines
 - dry: **ROPRON**
 - wet: **SAFIR**



Dispersing Systems

principle:

dispersing medium	in measuring zone of sensor	
	turbulent	statical

dry	RODOS	GRADIS
size ranges	$1 < x < 2500 \mu\text{m}$	$100 < x < 3500 \mu\text{m}$
quantity	5.0 g - 100 g	5.0 g - 100 g

wet	SUCCELL	CUVETTE
size ranges	$0.1 < x < 875 \mu\text{m}$	$0.5 < x < 1000 \mu\text{m}$
quantity	0.5 g - 2.0 g	0.001 g - 0.5 g

Laser Sensors Basic Equipment

★ central unit

in BA-, KA-, LA-, VARIO-cabinet of stainless steel with super pressure ventilation, and

- optical bench system
- HE-NE-laser or laser diode
- 31 channel multielement detector incl. auto alignment system (autofocus)
- data acquisition unit MECB-HEL
- high speed fibre optic communication interface


★ system computer

- HP VECTRA 486/XX
- HELOS software package
- high speed fibre optic communication interface
- XXX MByte hard-disk
- 3.5" floppy-disk 1.44 MByte
- HP 14" super VGA colour VDU
- keyboard
- MS-DOS license

Advantages of Dry Dispersion with RODOS

- ★ **complete and proper dispersion of all kinds of dry powders to $< 0.5 \mu\text{m}$**
 - ⇒ sensitive to changes in the fine particle range
- ★ **adjustable dispersing forces**
 - ⇒ best adaptation to dispersibility of different products
- ★ **variable sample size**
 - ⇒ minimising the statistical sampling error
- ★ **high measurement frequency due to short cycle time**
 - ⇒ no cleaning, no rinsing, of sample container
- ★ **no chemical interaction between particles and dispersing fluid**
 - ⇒ no dissolution, no surface reaction, ...
- ★ **no reprocessing of contaminated liquid**
- ★ **guaranteed long life time of disperser**
 - ⇒ low operating cost
 - ⇒ maximum up-time

Dry Dispersing System RODOS



A = 783 mm
B = 364 mm
C = 240 mm

auto-line

cycle time

steps	auto-RODOS	auto-SUCELL
1. information in	5"	5"
2. initialisation	10"	10"
3. fill in	15"	60"
4. sampling & dispersing	10" - 30"	10" - 30"
5. take out	15" - 35"	60" - 120"
6. evaluation	5"	5"
7. information out	5" - 60"	5" - 60"
<hr/>		
Σ 1. - 7.: [minutes]	< 2	5
<hr/>		
8. stand-by	1" - ∞	1" - ∞


auto-RODOS module

automatic dry dispersing system for use with HELOS laser diffraction analyser

- ★ **plug-in module for direct installation in HELOS measuring zone**
- ★ **integrates all functions for automatic particle size analysis, e.g. in cooperation with robots**
 - ☆ **mechanics:**
 - fully hardened dispersing line
 - polished inlet-funnel
 - ☆ **electronics:**
 - separate box on backside, fail-safe design
 - integrated security functions
 - local/remote operation
 - ☆ **integrated ancillary services:**
 - cover to eliminate external light, incl. operational panel
 - filters for compressed air, incl. monitoring of primary pressure
 - lens cleaning system with specially cleaned compressed air
 - streamline extraction unit, incl. monitoring of vacuum
 - fixed position relative to laser beam

Dry Dispersing System

auto-RODOS module






HELOS / MAGIC

laser diffraction system with automatic change of the measuring range (auto-ranging), including auto-focus features:


- ★ lense rotator for up to 4 lenses/ranges
- ★ detector positioning unit
- ★ MAGIC electronic
- ★ MAGIC software

⇒ MAGIC module cooperates with all types of auto-line dispersing units






21




22



Long Time Stability

Example auto-RODOS:



Blaine-values
$$S_M = \frac{6f}{\rho} \sum_i \frac{q_3(x_i)}{x_i} \Delta x_i$$
 (fine tail is significant)

- ★ standard deviation: < 1%, typical 0,4%
- ★ Blaine values: $\pm 40 \text{ cm}^2/\text{g}$ for $4000 \text{ cm}^2/\text{g}$
- ★ life time: ca. 250kg PZ35F
i. e. 25000 samples, 10 g each or typical 1 year