# Measuring sprays

## Particle size analysis: a routine procedure

*Because in spray products the active ingredient is delivered in the form of droplets the size of those droplets is critical for the product's efficacy and application.*

When developing aerosol and spray products for cosmetics there are two factors to be considered when deciding on particle size:

1. What is the ideal particle size for the spray's specific purpose? For example hairsprays require a finer distribution, deodorants can be somewhat coarser.

2. If the spray is not to be ingested into the body then droplets between about 1 and 10 microns in size must be avoided, in line with the Johannesburg Convention, to ensure that they do not reach the lungs.

Particle size in sprays is dependent on several factors. The formulation plays a major part, as does the geometry of the nozzle or actuator. The pressure in the can also affects the


particle size of sprays using propellant gases: in pump sprays the way that the spray is triggered has an effect.

### Measurement of droplets

Particle size analysis using laser diffraction provides the ability to measure droplet size with a high degree of accuracy. The following three aspects should be observed when this is being done:

→ Product dispersion must be carried out in a reproducible way. Only readings taken under the same, precisely defined conditions, can be compared. This is of particular relevance when firing pump sprays.

▼ Measuring system for the particle size

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→ The area in which the measurement is taking place must be of a certain minimum width in order that a normal spray actuation can be measured and also that droplets do not impair the performance of the optics in the laser device.

→ Analysis of the readings should not be parameter-dependent, and should give data on the number of particles in different size-categories, as well as the being time related, in order to give a detailed dynamic picture of the of the particle size, number and position over time.

The instrument manufacturers, Sympatec of Clausthal-Zellerfeld, have developed a dispersion module for their Helos laser diffraction sensor. The module, called the Spray Sizer, holds the spray package, which is operated by a robotic firing mechanism. It presses down with a precisely adjustable force, stroke length and acceleration – all factors which have a ma-

ior effect on the particle size produced by spray pumps. These actuation parameters are fed into the windows-based operator software and converted by the device into operation settings, with no physical intervention by the user. It is only in this way that exactly reproducible results can be achieved.

The Helos laser diffraction sensor can measure particles with the highest degree of resolution and accuracy, with almost all of the droplets passing in front of the broad laser beam and the evaluation, carried out in line with Fraunhofer recommendations, can be made without the need to apply compensation factors. Temporal resolution of 2000 individual readings per second, i.e once every 0.5 milliseconds, allows a detailed picture of the spray shot to be drawn and the best settings for the spray actuation to be determined.

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