

(side panel removed)

exhaust

Ø 32

25

240

compressed air supply

power supply and communication

128

(rear panel removed)

168

128

1017

134

258

product inlet

159


133

326

905

Maßstab	Abt		Bearbeitet	Geprüft		SYMPA System-Partikel TC Technik G.m.b.H.
1:10	501	Datum	18.01.07			
Benennung						Zeichnungs-Nr.
MYTOS & VIBRI Modul						A525 750/2

Copyright © 2007 by Sympatec GmbH, System-Partikel-Technik, Am Pulverhaus 1, D-38678 Clausthal-Zellerfeld
 Weitergabe sowie Vervielfältigung dieser Unterlage, Verwertung und Mitteilung ihres Inhaltes nicht gestattet, soweit nicht ausdrücklich zugestanden. Zuwiderhandlungen verpflichten zu Schadenersatz. Alle Rechte für den Fall der Patenterteilung oder Gebrauchsmuster Eintragung vorbehalten.
 Copying of this document, and giving it to others and the use or communication of the contents there, are forbidden without express authority. Offenders are liable to the payment of damages. All rights are reserved of the grant of a patent or the registration of a utility model of design.



Sympatec GmbH
 System-Partikel-Technik
 Clausthal-Zellerfeld,
 den 18.01.07